placeholder text

UC1: Start System		
Goal	To initiate the system and prepare to scan	
Initiation	Operator opens the software on the computer.	
Precondition	The computer has been booted	
Postcondition	The software has activated Robot arm, 3D camera and Ultra-	
	sound scanner.	
Main Scenario	1. Operator starts System	
Exceptions		

UC2: Detect chest area		
Goal	To construct a depth image of Patients chest area, so that UC3	
	can be initiated	
Initiation	Operator selects 3D Scan	
Precondition	UC1 has been completed	
Postcondition	Patients chest area has been scanned	
Main Scenario		
	1. 3D camera scans chest area	
	2. 3D camera feeds boundary points to System	
	3. Operator verifies depth image	
	• [Depth image is distorted]	
Exceptions		
	• [Depth image is distorted]	
	1. Enter UC2, item 1.	

UC3: Detect chest area		
Goal	To perform a medical sonography of the chest area	
Initiation	Operator selects Ultrasound Scan	
Precondition	UC2 has been completed	
Postcondition	System has instructed Robot arm to move Ultrasound scanner	
	in order to construct a medical sonographic image of Patients	
	chest area	
Main Scenario		
	1. Operator starts Ultrasound scanner	
	2. Operator chooses 'Ultrasound Scan' in System	
	3. System utilizes depth image from UC2 to feed Robot arm positions and rotations in order to move Ultrasound scanner	
	4. System stops Robot arm	
	• [Operator stops Robot arm]	
	5. Operator stops Ultrasound scanner	
Exceptions		
	• [Operator stops Robot arm]	
	1. Operator chooses 'Stop Ultrasound Scan' in System.	
	2. System stops moving Robot arm prematurely.	

UC4: Stop system		
Goal	To stop System	
Initiation	Operator	
Precondition	UC1 has been completed	
Postcondition	System has shut down	
Main Scenario	<ol> <li>Operator closes System software by hitting the Close button in the top right corner</li> <li>Upon exit, Robot arm is set to default position</li> <li>System shuts down</li> </ol>	
Exceptions		